

Last Updated: October 22nd, 2019

AIT Laboratories
a  HealthTRACKRx company

1500 Interstate 35W, Denton, TX 76207
P: (940) 591 8444 | www.healthtrackrx.com

Table of Contents

AIT Laboratories HL7 Interface Standard Specifications v.2 (v.2.3)	3
1. Purpose	3
2. Overview	3
2.1 Statement of Intent	3
2.2 Scope	3
3. General Specifications	3
3.1 Communication Protocol	3
3.2 Application Processing Procedures	3
3.3 HL7 messages supported by CLIN1 LIS	4
3.4 Segments	4
3.5 Segment Data Fields	4
3.5.1 Order Message (ORM)	5
3.5.2 Result Message (ORU)	10
3.6.VALUE TYPE	Error! Bookmark not defined.
4. Communication Requirements	Error! Bookmark not defined.

AIT Laboratories HL7 Interface Standard Specifications v.2 (v.2.3)

I. Purpose

This document specifies an interface to CLIN1 LIS, HealthTrackRx Laboratory Information System (**LIS**) based upon the Health Level Seven (HL7) Standard. This interface forms the basis for the exchange of healthcare information between the CLIN1 LIS with all other Electronic Medical Record systems (**EMR**) and Laboratory Information systems (**LIS**).

2. Overview

2.1 Statement of Intent

HealthTrackRx developed and implemented a standardized interface in accordance with the HL7 Standard, used by the CLIN1 software application in communicating with other EMR and LIS to exchange healthcare information. The interface strictly adheres to the HL7 Standard. This interface specification is subject to modification and revision to incorporate changes, improvements, and enhancements. This version supports HL7 V 2.3 Standard.

2.2 Scope

This document describes messages transmitted between the CLIN1 LIS application via the interface engine, MIRTH (NEXTGEN) with the client's EMR or LIS. The purpose of these messages is to exchange healthcare information concerning laboratory orders and results for further medical services.

3. General Specifications

3.1 Communication Protocol

The HL7 Standard defines only the seventh level of the Open System Inter-connect (OSI) Protocol. This is the application level. Levels 1-6 involve primarily communication protocols. Refer to the HL7 Standard for some guidance in this area.

3.2 Application Processing Procedures

The HL7 Standard describes the basic procedures for application processing by the sending and receiving systems. Anyone wishing to interface with CLIN1 LIS should be familiar with the HL7 Standard V. 2.3

3.3 HL7 messages supported by CLIN1 LIS

ORM	Observation Orders	Incoming
ORU	Observation Results	Outgoing

3.4 Segments

The following HL7 segments are used to support the exchange of Laboratory information.

Segments	Definitions
MSH	Message Header
PID	Patient Identification
IN1	Insurance
GT1	Guarantor
DG1	Diagnosis
ORC	Common Order
OBR	Observation Request
OBX	Observation
NTE	Notes and Comment

3.5 Segment Data Fields

The segment definition tables list and describe the data fields in the segment and characteristics of their usage. The following information is specified about each data field.

Fields	Description
Sequence Number (SEQ)	The ordinal position of the data field within the segment. This number is used to refer to the data field in the text comments that follow the segment definition table.
Length (LEN)	The maximum numbers of characters that one occurrence of the data field may occupy.
Data Type (DT)	Restrictions on the contents of the data field as defined by the HL7 Standard.
Optionally (R/O/C)	Whether the data field is required, optional, or conditional in a segment. The designations are R - required; O (null) -optional; and C - conditional on the trigger event.

Repetition (RP/#):	Whether the field may repeat. The designations are N (null) - for no repetition allowed; Y - the field may repeat an indefinite or site determined number of times; and (integer) - the field may repeat up to the number of times specified in the integer.
Element Name	Globally unique descriptive name for the field

3.5.1 Order Message (ORM)

I. Segment MSH-Message Header

This segment defines the intent, source, destination and some specifics of the syntax of a message

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	1		R	FIELD SEPARATOR	
2	4		R	ENCODING CHARACTERS	
3	15		R	SENDING APPLICATION	
4	20		R	SENDING FACILITY	ACCOUNT NUMBER
5	30		R	RECEIVING APPLICATION	
6	30		R	RECEIVING FACILITY	
7			R	DATE/TIME OF MESSAGE	
9	7		R	MESSAGE TYPE	
10			R	DATE/TIME OF MESSAGE	
11	1		R	PROCESSING ID	P-PRODUCTION; T-TESTING
12	8		R	VERSION ID	v.2.3

II. Segment: PID- Patient Identification

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	2	ID	R	SEQUENCE NUMBER	
2	15	NM	R	PATIENT ID (External ID)	

3	15	NM	R	PATIENT ID (Internal ID)	AIT patient Identifier
5	48	ST	R	PATIENT NAME	5.1-LastName; 5.2-FirstName
7		DT	R	DATE OF BIRTH	
8	1	ST	R	SEX	
11	106	XAD	R	ADDRESS	11.1-Address1; 11.2-Address2; 11.3-City; 11.4-State, 11.5-ZipCode

III. Segment: GT1- Guarantor

GT1 segment is required for Third Party Billing.

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	4	SI	R	SET ID	
2	250	CX	O	GUARANTOR NUMBER	
3	48	ST	R	GUARANTOR NAME	3.1-LastName; 3.2-FirstName
4	48	ST	O	GUARANTOR SPOUSE NAME	4.1-LastName; 4.2-FirstName
5	106	XAD	R	GUARANTOR ADDRESS	5.1-Address1; 5.2-Address2; 5.3-City; 5.4-State, 5.5-ZipCode
6	40	ST	O	GUARANTOR PHONE NUMBER	
8	28	DT	R	GUARANTOR DOB	
9	1	ST	R	GUARANTOR SEX	
10	2	ST	O	GUARANTOR TYPE	
11	80	CE	R	GUARANTOR RELATIONSHIP	
12	11	ST		GUARANTOR SSN	

IV. Segment: PV1 – Patient Visit

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	4	SI	R	SET ID	
2	250	CX	O	PATIENT CLASS	

3	48	ST	O	ASSIGNED PATIENT LOCATION	
4	48	ST	O	ADMISSION TYPE	
5	106	XAD	O	PREADMIT NUMBER	
6	40	ST	O	PRIOR PATIENT LOCATION	
7			R	ORDERING PROVIDER	7.1-Provider's NPI; 7.2-Provider's LastName; 7.3-Provider's FirstName
8	28	DT	O	REFERRING DOCTOR	
9	1	ST	O	CONSULTING DOCTOR	

V. Segment: IN1- Insurance

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	4	SI	R	SET ID	
2	60	CE	O	INSURANCE PLAN ID	
3	60	CX	R	INSURANCE COMPANY ID	
4	130	ST	R	INSURANCE COMPANY NAME	
5	106	XAD	O	INSURANCE COMPANY ADDRESS	5.1-Address1; 5.2-Address2; 5.3-City; 5.4-State, 5.5-ZipCode
6	48	ST	O	INSURANCE COMPANY CONTACT PERSON	6.1-LastName; 6.2-FirstName
7	40	ST	O	INSURANCE COMPANY PHONE NUMBER	
8	12	XON	O	GROUP NUMBER	
16	48	ST	R	NAME OF INSURED	8.1-LastName; 8.2-FirstName
17	80	CE	R	INSURED RELATIONSHIP TO PATIENT	
18	26	DT	R	INSURED DOB	
19	160	XAD		INSURED ADDRESS	19.1-Address1; 19.2-Address2; 19.3-City; 19.4-State, 19.5-ZipCode
22	2	ST	R	COORDINATION OF BENEFITS PRIORITY	

36	15	ST	R	POLICY NUMBER	
43	1	IS	R	INSURED SEX	

VI. Segment DG1 – Diagnosis

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	4	ID		SET ID	
2	2	ID	R	DIAGNOSIS CODING METHOD	
3	250	CE	R	DIAGNOSIS CODE	
4	40	ST	R	DIAGNOSIS DESCRIPTION	

*Required on all orders

VII. Segment ORC- Common Order

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	2	ID		SEQUENCE NUMBER	
2	40	EI	R	PLACER ORDER NUMBER	Receiving EMR/Lab Accession Number
3	10	EI	R	FILLER ORDER NUMBER	AIT Lab Accession Number
9	26	TS	R	DATE/TIME OF TRANSACTION	
12	50	XCN	R	ORDERING PROVIDER	16.1-Provider's NPI; 16.2-Provider's LastName; 16.3-Provider's FirstName

VIII. Segment: OBR -Observation Request

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	2	SI		SEQUENCE NUMBER	
2	40	EI	R	PLACER ORDER NUMBER	Receiving Lab Accession Number
3	10	EI	R	FILLER ORDER NUMBER	AIT Lab Accession Number
4	40	CE	R	OBSERVATION IDENTIFIER	4.1-Procedure Code; 4.2-Prodecue Code Name

6		TS	R	REQUESTED DATE/TIME	Order DateTime
7		TS	R	OBSERVATION DATE/TIME	Order DateTime
14		TS		SPECIMEN RECEIVED DATE/TIME	
15		SPS		SPECIMEN SOURCE	
16	50	XCN	R	ORDERING PROVIDER	16.1-Provider's NPI; 16.2-Provider's LastName; 16.3-Provider's FirstName

IX. Segment OBX- Observation (AOE only)

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	2	SI		SEQUENCE NUMBER	
2	2	ID		VALUE TYPE	Refer to 3.6 Value Type
3	40	CE	R	OBSERVATION IDENTIFIER	3.1-Procedure Code; 3.2-Procedure Name
4		VARIES	R	OBSERVATION SEQUENCE NUMBER	
5	250	CE	R	OBSERVATION VALUE	

X. Segment NTE- Notes

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	4	SI	O	SEQUENCE NUMBER	
2	8	ID	O	SOURCE OF COMMENT	
3	MAX	PT	O	COMMENT	

XI. Sample Orders

Orders 1

MSH|^~&||CLIN1|HTRX|20190606111159||ORM^O01|10002|P|2.3||NE|NE

PID|1||1812|HTRXTEST^TEST||19991108|M||1010 Easy St^^Anaheim^CA^13123|(555) 555-1212||||19-022570|

PV1|1||^&||||11000000000^DRTEST^TEST^|||||||4341|||||||20190511145900|20190606111159|

GT1|1||HTRXTEST^TEST||1010 Easy St^^Anaheim^CA^13123|(555) 555-1212||19991108||01|||||||

IN1|1||87726|UnitedHealth Care (WAC)|PO Box 30755^^Salt Lake City^UT^841300755||877-842-3210|9152587|||20190301|||HTRXTEST^TEST|01|19991108|1010 Easy St^^Anaheim^CA^13123||1|||||||||||||||||T

ORC|NW|0000003009|||20190606070000|11000000000^DRTEST^TEST^|11000000000^DRTEST^TEST^^^^^U|^^^&|

OBR|1|0000003009||127UALLZ^COMPREHENSIVE SCREEN AND REFLEX PANEL||20190606070000||^|11000000000^DRTEST^TEST^|1^S|

OBX|1|ST|MEDLIST^MEDLIST||Atenolol, Lisinopril, Naproxen, Cyclobenzaprine PRN, Mirtazapine (Remeron) PRN, Hydroxyzine (Vistaril) PRN, Hydroxyzine (Vistaril) PRN,

DG1|1|I10|F10.20|Alcohol use disorder, severe|

DG1|2|I10|F41.1|Generalized Anxiety Disorder|

3.5.2 Result Message (ORU)

I. Segment MSH-Message Header*

This segment defines the intent, source, destination and some specifics of the syntax of a message

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	1		R	FIELD SEPARATOR	
2	4		R	ENCODING CHARACTERS	
3	15		R	SENDING APPLICATION	"MIRTH"
4	20		R	SENDING FACILITY	"HTRX"
5	30		R	RECEIVING APPLICATION	
6	30		R	RECEIVING FACILITY	
7			R	DATE/TIME OF MESSAGE	
9	7		R	MESSAGE TYPE	
10			R	DATE/TIME OF MESSAGE	
11	1		R	PROCESSING ID	P-PRODUCTION; T-TESTING
12	8		R	VERSION ID	v.2.3

II. Segment: PID- Patient Identification*

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	2	ID	R	SEQUENCE NUMBER	
2	15	NM	R	PATIENT ID (External ID)	
3	15	NM	R	PATIENT ID (Internal ID)	AIT patient Identifier

5	48	ST	R	PATIENT NAME	5.1-LastName; 5.2-FirstName
7		DT	R	DATE OF BIRTH	
8	1	ST	R	SEX	
11	106	XAD	R	ADDRESS	11.1-Address1; 11.2-Address2; 11.3-City; 11.4-State, 11.5-ZipCode

III. Segment ORC- Common Order*

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	2	ID		SEQUENCE NUMBER	
2	40	EI	R	PLACER ORDER NUMBER	Receiving EMR/Lab Accession Number
3	10	EI	R	FILLER ORDER NUMBER	AIT Lab Accession Number
9	26	TS	R	DATE/TIME OF TRANSACTION	
12	50	XCN	R	ORDERING PROVIDER	16.1-Provider's NPI; 16.2-Provider's LastName; 16.3-Provider's FirstName

IV. Segment: OBR -Observation Request*

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	2	SI		SEQUENCE NUMBER	
2	40	EI	R	PLACER ORDER NUMBER	Receiving Lab Accession Number
3	10	EI	R	FILLER ORDER NUMBER	AIT Lab Accession Number
4	40	CE	R	OBSERVATION IDENTIFIER	4.1-Procedure Code; 4.2-Procedure Code Name
6		TS	R	REQUESTED DATE/TIME	Order DateTime
7		TS	R	OBSERVATION DATE/TIME	Order DateTime
14		TS		SPECIMEN RECEIVED DATE/TIME	
15		SPS		SPECIMEN SOURCE	

16	50	XCN	R	ORDERING PROVIDER	16.1-Provider's NPI; 16.2-Provider's LastName; 16.3-Provider's FirstName
21	30	ST		PERFORMING LABORATORY LOCATION	22.1- Laboratory Name; 22.2- Address1; 22.3-Address2; 22.4-City; 22.5-State; 22.6-ZipCode
22		TS	R	REPORTING DATE/TIME	Date/Time at which the results were reported
25	2	ID	R	OBSERVATION RESULT STATUS	F-FINAL, C-CORRECTED
32	50			PRINCIPAL RESULT INTERPRETER	FULL NAME
34				TECHNICIAN	INITIALS

V. Segment OBX- Observation

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	2	SI		SEQUENCE NUMBER	
2	2	ID		VALUE TYPE	Refer to 3.6 Value Type
3	40	CE	R	OBSERVATION IDENTIFIER	3.1-Procedure Code; 3.2-Procedure Name
5		VARIES	R	OBSERVATION VALUE	For embedded PDF: OBX.5.3 - 22639-9^PDF Report OBX.5.4 - ^^^^Base64^Embedded PDF Data
6	250	CE	R	UNITS	
7	60	ST		REFERENCE RANGES	
8	250	IS	R	ABNORMAL FLAGS	
11	1	ID	R	OBSERVATION RESULT STATUS	F-FINAL, C-CORRECTED
14		TS	R	RECEIVING DATE/TIME	
15	250	CE		PRODUCER'S ID	Performing Location ID
16	50	XCN		RESPONSIBLE OBSERVER	Laboratory Analyst Initials
23	50	ST		PERFORMING LABORATORY NAME	
24	50	ST		PERFORMING LABORATORY LOCATION	22.1- Address1; 22.2-Address2; 22.3-City; 22.4-State; 22.5-ZipCode

VI. Segment NTE- Notes

SEQ	LEN	DT	R/O/C	ELEMENT NAME	COMMENTS
1	4	SI	O	SEQUENCE NUMBER	
2	8	ID	O	SOURCE OF COMMENT	
3	MAX	PT	O	COMMENT	

VII. Sample Results

Results 1

MSH|^~\&|CLIN1LIS|HTRX|||20190910160823||ORU^R01|20190910160823|P|2.3|||AL|AL

PID|1|005590854|0000469252||HTRXTEST^CINDY^||19981214|F|||||||||

ORC|RE|5666333^|10956083||CM|||||20190910160823|||1000000000^DRTEST^DANIEL|||||||

OBR|1|5666333|10956083|CHNUC^URINE

CONFIRM^L|Routine|201909101530|201909101530|||||201910101556|URINE|1000000000^DRTEST^DANIEL
||CHNUC|||HTRX^1500 interstate 35 West^^Denton^TX^76207^|20190910000000|||F|||||||PUPL|

OBX|1|ST|MEDLIST^MEDLIST||ALPRAZOLAM|||||F|||||||||

NTE|1||ALPRAZOLAM LAST DOSE TAKEN 10/10/19 3:23 PM||

OBX|2|NM|LORAZ^Lorazepam||0.0|ng/mL|0.1-

50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|3|ST|LORAZ^Lorazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500
interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|4|NM|DIAZP^Diazepam||0.0|ng/mL|0.1-50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500
interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|5|ST|DIAZP^Diazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500
interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|6|NM|NORDIAZ^Nordiazepam||0.0|ng/mL|0.1-

50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|7|ST|NORDIAZ^Nordiazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|
1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|8|NM|OXAZEPAM^Oxazepam||0.0|ng/mL|0.1-

50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|9|ST|OXAZEPAM^Oxazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1
500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|10|NM|TEMAZ^Temazepam||0.0|ng/mL|0.1-
50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|11|ST|TEMAZ^Temazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|15
00 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|12|NM|CLONAZP^Clonazepam||0.0|ng/mL|0.1-
50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|13|ST|CLONAZP^Clonazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|
1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|14|NM|CLONAZ^Clonazepam metabolite||0.0|ng/mL|0.1-
50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|15|ST|CLONAZ^Clonazepam
metabolite||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|16|NM|FLUNIP^Flunitrazepam||0.0|ng/mL|0.1-
50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|17|ST|FLUNIP^Flunitrazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|
1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|18|NM|FLUNI^Flunitrazepam metabolite||0.0|ng/mL|0.1-
50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|19|ST|FLUNI^Flunitrazepam
metabolite||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|20|NM|ALPRAZP^Alprazolam||34.0|ng/mL|0.1-
50.0|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|21|ST|ALPRAZP^Alprazolam||Below Reporting
Limit.|ng/mL|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|22|NM|ALPRAZ^Alprazolam metabolite||122.0|ng/mL|0.1-
50.0|H|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|23|ST|ALPRAZ^Alprazolam metabolite||POSITIVE Consistent with
prescriptions.|ng/mL||H|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|24|NM|TRIAZP^Triazolam||0.0|ng/mL|0.1-
50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|25|ST|TRIAZP^Triazolam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|26|NM|TRIAZP^Triazolam metabolite||0.0|ng/mL|0.1-50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|27|ST|TRIAZP^Triazolam metabolite||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|28|NM|FLURAZP^Flurazepam||0.0|ng/mL|0.1-50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|29|ST|FLURAZP^Flurazepam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|30|NM|FLURAZP^Flurazepam metabolite||0.0|ng/mL|0.1-50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|31|ST|FLURAZP^Flurazepam metabolite||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|32|NM|CHLRDIAZ^Chlordiazepoxide||0.0|ng/mL|0.1-50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|33|ST|CHLRDIAZ^Chlordiazepoxide||NEGATIVE|ng/mL||L||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|34|NM|CLOBAZ^Clobazam||0.0|ng/mL|0.1-50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|35|ST|CLOBAZ^Clobazam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|36|NM|ESTAZ^Estazolam||0.0|ng/mL|0.1-50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|37|ST|ESTAZ^Estazolam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|38|NM|MIDAZP^Midazolam||0.0|ng/mL|0.1-50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|39|ST|MIDAZP^Midazolam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|40|NM|MIDAZ^Midazolam metabolite||0.0|ng/mL|0.1-50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|41|ST|MIDAZ^Midazolam metabolite||NEGATIVE|ng/mL||L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

Results 2 (embedded PDF)

MSH|^~\&|CLIN1LIS|HTRX|||20190910160823||ORU^R01|20190910160823|P|2.3|||AL|AL

PID|1|005590854|0000469252||HTRXTEST^CINDY^||19981214|F|||||||||

ORC|RE|895562^|10956083||CM|||20190910160823||1000000000^DRTEST^DANIEL|||||||

OBR|1|895562|10956083|CHNUB^CHN URINE
CONFIRM^L|Routine|201909101530|201909101530|||||201910101556|URINE|1000000000^DRTEST^DANIEL
||CHNUB|||HTRX^1500 interstate 35 West^^Denton^TX^76207^|20190910000000||F|||||||PUPL|

OBX|1|ST|MEDLIST^MEDLIST||ALPRAZOLAM|||||F|||||||||

NTE|1|I|ALPRAZOLAM LAST DOSE TAKEN 10/10/19 3:23 PM||

OBX|2|NM|LORAZ^Lorazepam||0.0|ng/mL|0.1-50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|3|ST|LORAZ^Lorazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|4|NM|DIAZP^Diazepam||0.0|ng/mL|0.1-50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|5|ST|DIAZP^Diazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|6|NM|NORDIAZ^Nordiazepam||0.0|ng/mL|0.1-50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|7|ST|NORDIAZ^Nordiazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|8|NM|OXAZEPAM^Oxazepam||0.0|ng/mL|0.1-50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|9|ST|OXAZEPAM^Oxazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|10|NM|TEMAZ^Temazepam||0.0|ng/mL|0.1-50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|11|ST|TEMAZ^Temazepam||NEGATIVE|ng/mL||L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|12|NM|CLONAZP^Clonazepam||0.0|ng/mL|0.1-
50.0|L||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|13|ST|CLONAZP^Clonazepam||NEGATIVE|ng/mL||L||F|||201909101520|HTRX|pupadhyay|||||HTRX|
1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|14|NM|CLONAZ^Clonazepam metabolite||0.0|ng/mL|0.1-
50.0|L||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|15|ST|CLONAZ^Clonazepam
metabolite||NEGATIVE|ng/mL||L||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|16|NM|FLUNIP^Flunitrazepam||0.0|ng/mL|0.1-
50.0|L||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|17|ST|FLUNIP^Flunitrazepam||NEGATIVE|ng/mL||L||F|||201909101520|HTRX|pupadhyay|||||HTRX|
1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|18|NM|FLUNI^Flunitrazepam metabolite||0.0|ng/mL|0.1-
50.0|L||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|19|ST|FLUNI^Flunitrazepam
metabolite||NEGATIVE|ng/mL||L||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|20|NM|ALPRAZP^Alprazolam||34.0|ng/mL|0.1-
50.0|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|21|ST|ALPRAZP^Alprazolam||Below Reporting
Limit.|ng/mL|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|22|NM|ALPRAZ^Alprazolam metabolite||122.0|ng/mL|0.1-
50.0|H||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|23|ST|ALPRAZ^Alprazolam metabolite||POSITIVE Consistent with
prescriptions.|ng/mL||H||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|24|NM|TRIAZP^Triazolam||0.0|ng/mL|0.1-
50.0|L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|25|ST|TRIAZP^Triazolam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500
interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|26|NM|TRIAZ^Triazolam metabolite||0.0|ng/mL|0.1-
50.0|L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|27|ST|TRIAZ^Triazolam
metabolite||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|28|NM|FLURAZP^Flurazepam||0.0|ng/mL|0.1-
50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|29|ST|FLURAZP^Flurazepam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1
500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|30|NM|FLURAZ^Flurazepam metabolite||0.0|ng/mL|0.1-
50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|31|ST|FLURAZ^Flurazepam
metabolite||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|32|NM|CHLRDIAZ^Chlordiazepoxide||0.0|ng/mL|0.1-
50.0|L|||F|||201909101520|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|33|ST|CHLRDIAZ^Chlordiazepoxide||NEGATIVE|ng/mL||L||F|||201909101520|HTRX|pupadhyay|||||
HTRX|1500 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|34|NM|CLOBAZ^Clobazam||0.0|ng/mL|0.1-
50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|35|ST|CLOBAZ^Clobazam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|150
0 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|36|NM|ESTAZ^Estazolam||0.0|ng/mL|0.1-
50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|37|ST|ESTAZ^Estazolam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500
interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|38|NM|MIDAZP^Midazolam||0.0|ng/mL|0.1-
50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|39|ST|MIDAZP^Midazolam||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|15
00 interstate 35 West^^Denton^TX^76207^|45D2009077|

OBX|40|NM|MIDAZ^Midazolam metabolite||0.0|ng/mL|0.1-
50.0|L|||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|41|ST|MIDAZ^Midazolam
metabolite||NEGATIVE|ng/mL||L||F|||201910101521|HTRX|pupadhyay|||||HTRX|1500 interstate 35
West^^Denton^TX^76207^|45D2009077|

OBX|41|RP|22639-9^PDF

[illegible]